

spirit of excellence






Distinctive design for an effective clean




Designed to be backloaded through the working channel. (insert brushless end into the luer connector, and remove from the distal end of the scope).

All of our brushes are validated for use by channel diameter. Take the guess work out of selecting a brush for cleaning by using the charts to determine which brushes go with each of our flexible fiberoptic and sensor scopes.

Brushes for flexible scopes









Trays for flexible scopes

Flexible Endoscope Tray, Metal (85850302)

- Sturdy, silicone brackets safely hold the contents in place during packaging, sterilization, and storage
- Stainless steel construction for protection during transport
- 22.5"L x 10.5"W x 2.9"D



Uro Flexible Endoscope Tray, Plastic (8585035)

- Sturdy, silicone brackets safely hold scope in place during packaging, sterilization and storage
- 21.3"L x 10"W x 3"H



Contact your Richard Wolf Sales Representative or call (800) 323-WOLF (9653) for ordering information.

