

RIWO > D-URS

RICHARD
WOLF 

spirit of excellence

Feel the difference.



The First 3-channel disposable sensor ureterorenoscope:

- Separate backflow channel in the instrument sheath
- 9.0 Fr. outer diameter and atraumatic 6.0 Fr. distal tip
- 3.6 Fr. working channel for auxiliary instruments
- 1.65 Fr. working channel for laser fibers

Combination of semi-rigid sheath and flexible distal tip capable of 300° bidirectional deflection. Unique instrument design enables completely new handling. Surgery can also be performed comfortably while seated.

For more detailed information, contact
your local Sales Representative or
Richard Wolf Customer Service.

Customer Service
(800) 323-WOLF (9653), option 2.
info@richardwolfusa.com

Re-thinking Ureteroscopy — Best of Both Worlds.